

BARIUM (ATOMIC ABSORPTION, FURNACE TECHNIQUE)

1.0 SCOPE AND APPLICATION

1.1 See Section 1.0 of Method 7000.

2.0 SUMMARY OF METHOD

2.1 See Section 2.0 of Method 7000.

3.0 INTERFERENCES

3.1 See Section 3.0 of Method 7000.

3.2 Barium is known to form a barium carbide in the graphite furnace. This less volatile carbide can cause losses of sensitivity and memory effects.

3.3 The long residence time and the high concentration of the analyte in the optical path of the graphite furnace can lead to severe physical and chemical interferences. Furnace parameters must be optimized to minimize these effects.

3.4 Because of possible chemical interaction, nitrogen should not be used as a purge gas.

3.5 Halide acids should not be used.

4.0 APPARATUS AND MATERIALS

4.1 For basic apparatus, see Section 4.0 of Method 7000.

4.2 Instrument parameters (general):

4.2.1 Drying time and temp: 30 sec at 125°C.

4.2.2 Ashing time and temp: 30 sec at 1200°C.

4.2.3 Atomizing time and temp: 10 sec at 2800°C.

4.2.4 Purge gas: Argon (nitrogen should not be used).

4.2.5 Wavelength: 553.6 nm.

4.2.6 Background correction: Not required.

4.2.7 Other operating parameters should be set as specified by the particular instrument manufacturer.

NOTE: The above concentration values and instrument conditions are for a Perkin-Elmer HGA-2100, based on the use of a 20- $\mu$ L injection, continuous-flow purge gas, and nonpyrolytic graphite. Smaller size furnace devices or those employing faster rates of atomization can be operated using lower atomization temperatures for shorter time periods than the above-recommended settings.

## 5.0 REAGENTS

5.1 See Section 5.0 of Method 7000.

5.2 Preparation of standards

5.2.1 Stock solution - Dissolve 1.7787 g barium chloride ( $\text{BaCl}_2 \cdot 2\text{H}_2\text{O}$ , analytical reagent grade) in water and dilute to 1 liter. Alternatively, procure a certified standard from a supplier and verify by comparison with a second standard.

5.2.2 Prepare dilutions of the stock solution to be used as calibration standards at the time of analysis. The calibration standards should be prepared using the same type of acid and at the same concentrations as in the sample after processing (0.5% v/v  $\text{HNO}_3$ ).

## 6.0 SAMPLE COLLECTION, PRESERVATION, AND HANDLING

6.1 See Chapter Three, Step 3.1.3, Sample Handling and Preservation.

## 7.0 PROCEDURE

7.1 Sample Preparation - The procedures for preparation of the sample are given in Chapter Three, Step 3.2.

7.2 See Method 7000, Step 7.3, Furnace Technique.

## 8.0 QUALITY ASSURANCE

8.1 See Section 8.0 of Method 7000.

## 9.0 METHOD PERFORMANCE

9.1 Precision and accuracy data are not available at this time.

## 10.0 REFERENCES

1. Methods for Chemical Analysis of Water and Wastes; U.S. Environmental Protection Agency. Office of Research and Development. Environmental Monitoring and Support Laboratory. ORD Publication Offices of Center for Environmental Research Information: Cincinnati, OH, 1983; EPA-600/4-79-020.

METHOD 7081  
BARIUM (ATOMIC ABSORPTION, FURNACE TECHNIQUE)

